The Thermo Scientific Nicolet 380 FT-IR spectrometer is a reliable, easy-to-use system that helps you *identify, quantify and verify* your samples. The instrument, software and accessories are seamlessly integrated to provide you with answers and solutions.

Thermo Scientific Nicolet 380 FT-IR Spectrometer

Proven and reliable FT-IR spectroscopy











The Nicolet[™] 380 spectrometer provides high-quality results and greater productivity through the use of:

- Smart components
- Automatic system and accessory recognition
- Fully compatible with standard and Smart accessories
- Full-size sample compartment for any type of transmission experiment
- Diagnostics and optimization
- Snap-in experiments
- Multi-media tutorials

At the spectrometer's core is proven Thermo Scientific Enhanced Synchronization Protocol (E.S.P.) technology, which brings integration and intelligence to every part of the Nicolet 380 spectrometer system that affects the success of your measurements – from setup to report generation. The spectrometer fits into a footprint the size of a typical desktop computer. Despite its compact size, the instrument offers a full-sized sample compartment that will hold all the accessories you need. Sampling efficiency with the Nicolet 380 spectrometer is fully realized when combined with Smart Accessories, which snap in place to customize and optimize the spectrometer for specific analyses.

Truly Intelligent Sampling

Smart Accessories provide an unprecedented level of communication with the Nicolet 380 spectrometer, with features including:

- Permanently aligned optics to guarantee optimal performance and ensure precision from one measurement to the next
- Rugged, enclosed design to protect optics
- Snap-in installation and automatic recognition of accessories to load your experiments and optimize the spectrometer for the sampling technique



- Diagnostic checks to ensure optimal performance and check for "human errors"
- Multimedia tutorials to help users get the most out of the accessory
- Automatic purge within the accessory to provide quick equilibration times when an accessory is changed, and exceptional measurement stability over the long term
- Spectral guality checks monitor data collection to ensure consistent results

In addition to Smart Accessories, the full-size sample compartment accommodates most commercially available accessories, helping you protect your investment by using your existing accessories.

Precise, Reliable Analysis

Thermo Fisher Scientific ISO 9001 engineering and manufacturing procedures and quality control process ensure reliable results. All components are pinned-in-place which means they are alignment-free and user-replaceable. The Nicolet 380 design allows you to achieve a wavelength precision better than 0.01 cm⁻¹ without using risky mathematical "fudge factors" that alter data and are difficult to explain in an SOP.

E.S.P. technology monitors the performance of all key components to ensure data quality. If the system detects something unusual, either the problem is fixed automatically or precise, easily understood recommendations are provided.

The laser beam goes straight into the interferometer, which contains the only moving part of the spectrometer - the moving mirror - which is dynamically aligned to compensate for any mirror tilt, shear, or sag.

System Verification

A comprehensive spectrometer qualification package is available to verify system performance in compliance with the widely-accepted ASTM E1421 protocols used for compliance with ISO 9001. The ASTM method may be run on a routine basis to verify the system is operating within acceptable performance levels, based on NIST-traceable calibrated spectral standards.



Optics


- Baseplate: single-piece, precision-cast and machined baseplate incorporating pinnedin-placed components for extra stability and reliability; standard purge ports and low purge volume design
- Mirror Optics: proprietary diamond-turned, pinned-in-place mirrors
- Layout: efficient design with short pathlength; a single mirror from the beamsplitter to the sample compartment; a single mirror from the sample compartment to the detector
- Source: long-lifetime, dual-mode, high-energy Ever-Glo source; pre-aligned, pinned source position for easy replacement
- Interferometer: rugged, frictionless electromagnetic drive; digital dynamic alignment and digital signal processing (DSP) control for long-term stability; AutoTune, our patented technology, for optimization of system throughput (U.S. Patent 5,883,712)
- Beamsplitter: proprietary Germanium on KBr (7800 - 350 cm⁻¹) beamsplitter optimized for high energy throughput in the infrared spectrum's "fingerprint" region
- Reference Laser: reference Helium Neon laser; pre-aligned, pinned-in-place, and user-replaceable
- Desiccant: reusable desiccant in easy-tochange replaceable bag; visible desiccant moisture indicator without opening the system cover
- Detector: high-performance DTGS detector mounted on pre-aligned, pinned-in-place baseplate for optimal sensitivity

Scientific Nicolet 380 Spectrometer Th

Accessories	Smart Accessories support	Standard
	Traditional accessory support	Standard
Performance	Spectral range	7800 – 350 cm ⁻¹ using proprietary KBr beamsplitter
	Optical resolution (apodized)	< 0.9 cm ⁻¹ resolution (standard)
	Peak-to-peak noise (1 minute scan)	< 2.2 x 10 ⁻⁵ Abs. (> 22,000:1)
	RMS noise (1 minute scan)	< 5.5 x 10 ⁻⁶ Abs.
	Wavenumber precision	Better than 0.01 cm ⁻¹ precision at 2000 cm ⁻¹
	Ordinate linearity (ASTM E1421)	< 0.1%T deviation from 0.0%T at 4 cm ⁻¹ resolution
Electronics	A/D converter	24-bit
	On-bench controller	DSP-based
Computer Interface	Operating system	Win 2000, Win XP, Vista
	Interface	USB 2.0
Warranty	1 year on complete system	
Spectrometer Dimensions	50 cm (w) x 58 cm (d) x 23 cm (h) with Smart Accessories; 29 cm (h) with full sample compartment	
Spectrometer Weight	24 kg	
Sample Compartment	Standard, full-sized; 21 cm (w) x 26 cm (d) x 20 cm (h)	
Regulatory Approvals		

©2004-2009 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.





PS50673 E 07/09M



Europe-Other +43 1 333 50 34 0 Italy +39 02 950 591 Finland/Norway/Sweden +46 8 556 468 00 **France** +33 1 60 92 48 00 Germany +49 6103 408 1014

India +91 22 6742 9434 Japan +81 45 453 9100 Latin America +1 608 276 5659 UK +44 1442 233555 Middle East +43 1 333 50 34 0 Netherlands +31 76 579 55 55

South Africa +27 11 570 1840 Spain +34 914 845 965 Switzerland +41 61 716 77 00 USA +1 800 532 4752 www.thermo.com



System Software

OMNIC Lite Software: menu and toolbar operation; fully customizable interface with password protection; integrated control of setup, collection, data manipulation, and reporting; comprehensive system diagnostics for the entire system and Smart Accessories; real-time display of spectral data and interactive parameter setting; context-sensitive help; multimedia tutorials for FT-IR theory, sampling, accessories, data collection, diagnostics, and other support

Software Options

OMNIC Lite software provides seamless interaction with optional add-on packages for advanced analysis, including:

- Macros Basic to automate standardized tasks and SOPs to a simple push-button procedure
- **TO Analyst** software for guantitative and advanced qualitative analysis of spectral data sets, including algorithms such as Discriminant Analysis and PLS
- OMNIC Specta for unsurpassed identification of unknown compounds, including multi-component search and contaminant analysis